

DILUWEL XL[®]

LARGE VOLUME GRAVIMETRIC DILUTOR

Combine speed and comfort !



LABGIG

 **ALLIANCE
BIO EXPERTISE**
Connected Microbiology

LARGE VOLUME GRAVIMETRIC DILUTOR

Combine speed and comfort !

DILUWEL XL®

Dilution time : 36 seconds for 3750 g

New record !




Motorised arm

Avoid cross contamination

7" touchpad

Gloves friendly



Colored display tower

Adjustable height (from 34 to 50 cm)



Precision balance

From 8.00 g to 5000.00 g



Remote pump

Optimize your workflow



Ideal under biosafety cabinets !



Waterproof ambidextrous control panel

Industrial quality



A MULTIFUNCTION DILUTER :

Go for a mode and choose the best holders for your applications !

LANDSCAPE MODE CONFIGURATIONS :



Multi-Container
Bottles / Bowls / Flasks



Tilt
400 ml bags



Classic
4000 ml bags

PORTRAIT MODE CONFIGURATIONS :



Multi-Container
Bottles / Bowls / Flasks



Pull and Tilt
400 ml bags



Classic
4000 ml bags



Ergonomic Workflow
4000 ml bags

VERSATILITY

3750 g Your dilution in **36 seconds**

1250 g in **12 seconds**

250 g in **8 seconds**

100 g in **7 seconds**



New record!


ADAPTABILITY



Compatible with all automated media and diluent preparators




DILUWEL XL® ↔ **MEDIAWEL®**
The perfect combo for maximum speed



ERGONOMY

Designed to avoid musculoskeletal disorders



FOOD SAMPLE PREPARATION


Blenders



Diluters




Tubes Fillers Pumps




AUTOMATED MEDIA PREPARATION

Media Preparators



Petri Fillers



SAMPLE COLLECTION

Air Samplers



BOOK YOUR DEMO !

<https://www.alliance-bio-expertise.com/en/store-locator/>



FOLLOW THE BLOG

alliance-bio-expertise.com



Address : Vidmahe Enclave, No 135, Ground floor,
Kodigehalli - Thindlu Main Road, opp. BESCOM
Office, IEHCS Layout, Post, Vidyaranyapura,
Bengaluru, Karnataka 560097

Email : info@labgig.in