

2023

Shenhua Science Technology Co., Ltd.

Electrophoresis Power Supply

Introduction

The universal electrophoresis instruments are suitable for horizontal chambers, vertical chambers, and transfer chambers. The high-definition color LCD screen provides clear and intuitive parameters and data.




Features

- High quality ABS material
- Streamlined design
- Voltage Accuracy 0.1 V, Current Accuracy 0.1mA
- Colorful LCD
- Auto micro current holding function
- Protection for overvoltage, overcurrent, overload and etc.

Parameters

Model	SH-DDY10P	SH-DDY20P
Channels	Parallel connection for 4 channels	
Output Voltage	6-600V	6-300V
Accurate	0.1V	
Output Current	4-600mA	5-2000mA
Accurate	0.1mA	
Weight/size	1.5kg/280x240x120mm	



Address : Vidmahe Enclave, No 135, Ground floor,
Kodigehalli - Thindlu Main Road, opp. BESCO
Office, IEHCS Layout, Post, Vidyaranyapura,
Bengaluru, Karnataka 560097

Email : info@labgig.in

